




Sentri™ Intersurgical EcoLite™ End Tidal CO₂ Mask

Delivers supplementary oxygen with CO₂ monitoring



Oxygen and Aerosol Therapy • Variable Oxygen Concentration (Low Flow)



Quality, innovation and choice



www.intersurgical.com

Intersurgical EcoLite™ Mask

Monitoring ETCO₂ is an invaluable tool for reliable, early warning of respiratory compromise. By using the Senti™ Intersurgical EcoLite™ ETCO₂ mask, no matter where in the hospital setting IV sedation is delivered, both the patient and clinician can take comfort knowing that early detection of inadequate ventilation is made easier with this innovative mask.

The Senti Intersurgical EcoLite Mask is a medium concentration oxygen mask with the addition of a sampling port. This port permits sampling of exhaled carbon dioxide from both the mouth and nose in non-intubated patients during the administration of supplementary oxygen.

Capnography is vital during sedation

The increased use of conscious sedation has created a need for a device to monitor respiratory depression.

The difference between conscious sedation and general anesthesia is sometimes very small. It is possible during conscious sedation that intravenous sedatives and narcotics administered to alleviate apprehension can result in the loss of consciousness and respiratory obstruction.

The Senti Intersurgical EcoLite Mask is simple to use: attach one end of the ETCO₂ monitoring line to a sidestream capnometer and the other to the luer lock port on the mask.

Incurved nose seal

Conforms to different nose shapes, designed to prevent oxygen entering patient's eyes

No metal nose clip

MRI compatible

Elastic can be positioned under or over the ears

Below ear position eliminates trauma to top of ears

Soft face seals

Increased patient comfort

Multi-channel oxygen tube

Oxygen still flows even if tube is kinked



Luer lock port

Secure ETCO₂ monitoring line connection



A choice of "under chin" or "on chin" positions

Two integral chin seals provide a better fit on a wider range of patient face shapes.



Comfortable for the environment

Various groups have questioned the use of PVC in medical products and its impact on the environment. Intersurgical's aim is to reduce the environmental impact of its products and processes. The utilization of new materials and the latest manufacturing technology has resulted in the elimination of PVC from the mask shell resulting in a reduced environmental impact.



More about our Eco Range

Make an inquiry

us.intersurgical.com/products/oxygen-and-aerosol-therapy/senti-etco2-medium-concentration-masks

Code	Description	Tube length	Box Qty.
1142015	Senti Intersurgical EcoLite, adult, mask with CO ₂ monitoring line and tube	2.1m	30
1141015 NEW	Senti Intersurgical EcoLite, adult, mask with oxygen tube	2.1m	30

IS10.24 US • Issue 6 03.18