



i-gel[®] pediatric range

The supraglottic airway with a non-inflatable cuff



Airway Management ▪ Airway Devices



Quality, innovation and choice

Interact with us



www.intersurgicalinc.com



Airway management has evolved

Introducing the pediatric i-gel®: a revolutionary single use supraglottic airway from Intersurgical.

i-gel® and natural airway management

The i-gel® is a truly unique single use, latex and PVC free airway device, representing the culmination of years of extensive research and development. Based on the original i-gel® design inspired by the physiology of the perilaryngeal framework itself, everything about the pediatric i-gel® has been designed to work in perfect unison with a child's anatomy, and the specific anatomical, physiological and pathological differences that make pediatric anesthesia particularly challenging.

Pediatrics

i-gel® is available in four pediatric as well as three adult sizes, making it applicable for use with patients from 2kg to 90+kg. The pediatric i-gel® is indicated in securing and maintaining a patent airway in routine and emergency anesthetics for operations of fasted patients during spontaneous or intermittent positive pressure ventilation (IPPV) and during resuscitation of the unconscious patient, by personnel who are suitably trained and experienced in the use of airway management techniques and devices.

i-gel® mirrors the anatomy

The shape, softness and contours accurately mirror the perilaryngeal anatomy to create the perfect fit. This innovative concept means that no cuff inflation is required. The i-gel® works in harmony with the patient's anatomy so that compression and displacement trauma are significantly reduced or eliminated.

The non-inflatable cuff

i-gel® gets its name from the soft gel-like material from which it is made. It is the innovative application of this material that has enabled the development of a unique non-inflatable cuff. This key feature means insertion of i-gel® is easy, rapid and consistently reliable.

The simple, safe and rapid solution

i-gel® is incredibly easy to use. A proficient user can achieve insertion of the i-gel® in less than 5 seconds. With no inflatable cuff, i-gel® provides a safe and rapid airway management solution.

Accurate and natural positioning

The i-gel® accurately and naturally positions itself over the laryngeal framework, providing a reliable perilaryngeal seal without the need for an inflatable cuff.

Additional information available

Further information regarding i-gel® including bibliographies, a clinical study database and product videos is available on the i-gel® website www.intersurgicalinc.com/info/igel.



Evidence

1. Beringer RM, Kelly F, Cook TM, Nolan J, Hardy R, Simpson T, White MC. A cohort evaluation of the paediatric i-gel® airway during anaesthesia in 120 children. *Anaesthesia* 2011 Dec;66(12):1121-6
2. Beylacq L, Bordes M, Semjen F, Cros AM. The i-gel®, a single-use supraglottic airway device with a non-inflatable cuff and an esophageal vent: an observational study in children. *Acta Anaesthesiol Scand.* 2009 Mar;53(3):376-9
3. Kim MS, Oh JT, Min JY, Lee KH, Lee JR. A randomised comparison of the i-gel® and the Laryngeal Mask Airway Classic™ in infants. *Anaesthesia.* 2014 Apr;69(4):362-7
4. Das B, Mitra S, Jamil SN, Varshney RK. Comparison of three supraglottic devices in anesthetised paralyzed children undergoing elective surgery. *Saudi J Anaesth.* 2012 Jul;6(3):224-8
5. Mitra S, Das B, Jamil SN. Comparison of Size 2.5 i-gel® with Proseal LMA™ in Anaesthetised, Paralyzed Children Undergoing Elective Surgery. *N Am J Med Sci.* 2012 Oct;4(10):453-7
6. Goyal R, Shukla RN, Kumar G. Comparison of size 2 i-gel® supraglottic airway with LMA-Pro-Seal™ and LMA-Classical™ in spontaneously breathing children undergoing elective surgery. *Paediatr Anaesth.* 2012 April; 22(4):355-9
7. Lee JR, Kim MS, Kim JT, Byon HJ, Park YH, Kim HS, Kim CS. A randomised trial comparing the i-gel® with the LMA Classic (TM) in children. *Anaesthesia.* 2012 Jun;67(6):606-11
8. L Theiler, R Greif. Clinical evidence for the use of the i-gel® pediatric anesthesia. *The i-gel® supraglottic airway: Medical procedures, testing and technology.* Nova Science Publishers, 2013
9. Y Abukawa, K Hiroki, M Ozaki. Evaluation of the i-gel® airway in children. *The i-gel® supraglottic airway: Medical procedures, testing and technology.* Nova Science Publishers, 2013

Find out more



Features and benefits

The i-gel[®] has a host of features that provide significant benefits to the patient and the clinician.



15mm connector

Reliable connection to any standard catheter mount or connection

Proximal end of gastric channel

Clearly displayed product information

For quick easy reference. Includes confirmation of size and weight guidance



Gastric channel

The i-gel[®] incorporates a gastric channel (except size 1). It provides an early warning of regurgitation, allows for the passing of a nasogastric tube to empty stomach contents and facilitates venting

Integral bite block

Reduces the possibility of airway channel occlusion

Buccal cavity stabilizer

Aids insertion and eliminates the potential for rotation

Epiglottic rest

Reduces the possibility of epiglottic 'down folding' and airway obstruction

The non-inflatable cuff

Made from a unique soft gel-like material allowing ease of insertion and reduced trauma

Distal end of gastric channel



Correct taping technique

It is important that as soon as insertion has been completed, the i-gel[®] is held until and while the device is secured in place.

The i-gel[®] should be taped in place 'maxilla to maxilla' (upper jaw to upper jaw) as shown.



Innovative packaging

The pediatric i-gel® supraglottic airway is supplied in a fully recyclable cage pack. This unique packaging protects the i-gel® in transit and ensures that it maintains its anatomical shape. i-gel® is available in four pediatric sizes.



Code	Description	Color	Size	Weight	Box Qty.
8225000	i-gel®, supraglottic airway	○	2.5 Large pediatric	25–35kg (55-77lbs)	10
8202000	i-gel®, supraglottic airway	●	2 Small pediatric	10–25kg (22-55lbs)	10
8215000	i-gel®, supraglottic airway	●	1.5 Infant	5–12kg (11-25lbs)	10
8201000	i-gel®, supraglottic airway	●	1 Neonate	2–5kg (5-11 lbs)	10

[Make an inquiry](#)
[Find out more](#)
[Watch the video](#)

Sterile

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 224.2 / IS6.8 US • Issue 7 07.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan
Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA