




Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	QuadraLite masks
Neonatal	Preterm	< 2	0	 Code: 7191030
	Term	2 – 4		
Infant	3 months	6	1	 Code: 7191030
	6 months	8		
	9 months	9		
Pediatric	1 year	10	2	 Code: 7192030
	2 years	12		
	3 years	14		
	4 years	16		
	5 years	18		
Small Adult	6 years	20	3	 Code: 7193030
	7 years	23		
	8 years	26		
	9 years	28		
	10 years	32		
	11 years	35		
Medium Adult	12 years	39	4	 Code: 7193030
	13 years	44		
	14 years	50		
	15 years	54		
Large Adult	16 years	58	5	 Code: 7194030
	16+ years	60 – 120		
Extra Large Adult	16+ years	> 120	6	 Code: 7194030