

QuadraLite™

Four sizes. One great fit.
Manage your patient not the mask.



QuadraLite™

Soft seal

A comfortable airtight fit with flexibility for different face sizes

Transparent shell

For excellent visibility

Color coding

Color matching with other Intersurgical products for easy selection

Anatomical shape

Provides a precise fit combined with lower dead space

Less means more, to the power of four

The QuadraLite is a soft and comfortable anesthesia mask that is available in four sizes to fit most patients. The unique anatomical shape combined with its flexible seal means that only four sizes are required, unlike traditional cushion masks which often require six or seven sizes. This allows you to hold less stock and still provide a complete solution for most of your patients.

The greener, safer choice








QuadraLite masks contain no PVC and are phthalate free.



Comfort without the cushion

QuadraLite masks do not have a cushion, but an innovative use of two different materials which offer an optimum level of flexibility and comfort. A crystal-clear shell enhances patient visibility.

Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	QuadraLite masks
Neonatal	Preterm	< 2	0	 Code: 7191030
	Term	2 – 4		
Infant	3 months	6	1	 Code: 7191030
	6 months	8		
	9 months	9		
Pediatric	1 year	10	2	 Code: 7192030
	2 years	12		
	3 years	14		
	4 years	16		
	5 years	18		
Small Adult	6 years	20	3	 Code: 7193030
	7 years	23		
	8 years	26		
	9 years	28		
	10 years	32		
	11 years	35		
	12 years	39		
Medium Adult	13 years	44	4	 Code: 7193030
	14 years	50		
	15 years	54		
	16 years	58		
Large Adult	16+ years	60 – 120	5	 Code: 7194030
Extra Large Adult	16 + years	> 120	6	 Code: 7194030

Ordering information



<div><div>7194030</div><div>22ID</div></div> <div><div>7193030</div><div>22ID</div></div> <div><div>7192030</div><div>22ID</div></div> <div><div>7191030</div><div>15OD</div></div>			
QuadraLite anesthesia mask			
7194030	Large anesthetic face mask with hook ring size 4		50
7193030	Medium anesthetic face mask with hook ring size 3		50
7192030	Small anesthetic face mask with hook ring size 2		50
7191030	Infant anesthetic face mask size 1		50