



# Solus™ laryngeal mask airways

A comprehensive range of single use laryngeal mask airways



Airway Management ▪ Airway Devices



Quality, innovation and choice

Interact with us



[www.intersurgicalinc.com](http://www.intersurgicalinc.com)



## Solus™ laryngeal mask airways

A comprehensive range of single use laryngeal mask airways suitable for anesthesia and emergency medicine.

Designed to meet the needs of a variety of clinical scenarios, all Solus laryngeal mask airways are single use, latex free and supplied sterile.

### Low friction material

Atraumatic insertion & removal

### Classic cuff shape

Renowned for reliability

### Essential information

Clearly visible for quick reference




### Integral inflation line

Reduces the risk of damage whilst in use

### High quality valve

To ensure continual cuff integrity



[Watch the video](#)

### Solus Standard

Using specially chosen materials, the Solus Standard is manufactured with the quality, flexibility and strength required to work effectively in a range of demanding clinical applications, including pediatric anesthesia. Available in a full size range, from adult to neonate.



[Make an inquiry](#)

Code	Description	Size	Weight	Box Qty.
8005000	Solus Standard, laryngeal mask airway	5 Large adult	70+kg (154lbs)	20
8004000	Solus Standard, laryngeal mask airway	4 Medium adult	50–70kg (110–154lbs)	20
8003000	Solus Standard, laryngeal mask airway	3 Small adult	30–50kg (65–110lbs)	20
8025000	Solus Standard, laryngeal mask airway	2.5 Large pediatric	20–30kg (44–65lbs)	20
8002000	Solus Standard, laryngeal mask airway	2 Small pediatric	10–20kg (22–44lbs)	20
8015000	Solus Standard, laryngeal mask airway	1.5 Infant	5–10kg (11–22lbs)	20
8001000	Solus Standard, laryngeal mask airway	1 Neonate	<5kg (<11lbs)	20



## Solus Satin

Solus Satin is our range of single use laryngeal mask airways featuring a softer airway tube providing increased flexibility.



Make an inquiry








Code	Description	Size	Weight	Box Qty.
8005003	Solus Satin, laryngeal mask airway	5 Large adult	70+kg (154lbs)	20 
8004003	Solus Satin, laryngeal mask airway	4 Medium adult	50–70kg (110–154lbs)	20 
8003003	Solus Satin, laryngeal mask airway	3 Small adult	30–50kg (65–110lbs)	20 

## Solus MRI-safe

A range of single use, sterile, laryngeal mask airways fitted with specially designed non-ferrous valves, guaranteed not to interfere with the magnet in an MRI scanner. The plastic valve has been selected to ensure full reliability throughout the shelf life of each device.



Make an inquiry

Code	Description	Size	Weight	Box Qty.
8005002	Solus MRI-safe, laryngeal mask airway	5 Large adult	70+kg (154lbs)	20 
8004002	Solus MRI-safe, laryngeal mask airway	4 Medium adult	50–70kg (110–154lbs)	20 
8003002	Solus MRI-safe, laryngeal mask airway	3 Small adult	30–50kg (65–110lbs)	20 
8025002	Solus MRI-safe, laryngeal mask airway	2.5 Large pediatric	20–30kg (44–65lbs)	20 
8002002	Solus MRI-safe, laryngeal mask airway	2 Small pediatric	10–20kg (22–44lbs)	20 
8015002	Solus MRI-safe, laryngeal mask airway	1.5 Infant	5–10kg (11–22lbs)	20 
8001002	Solus MRI-safe, laryngeal mask airway	1 Neonate	<5kg (<11lbs)	20 





### Solus Curve

Solus Curve is our laryngeal mask for those who prefer the insertion characteristics of a device with a curved airway tube.




Make an inquiry

Code	Description	Size	Weight	Box Qty.
8005050	Solus Curve, laryngeal mask airway	5 Large adult	70+kg (154lbs)	20 
8004050	Solus Curve, laryngeal mask airway	4 Medium adult	50–70kg (110–154lbs)	20 
8003050	Solus Curve, laryngeal mask airway	3 Small adult	30–50kg (65– 110lbs)	20 

### Solus Flexible

Solus Flexible provides an ideal solution for shared airway procedures such as ENT, dental, oro-maxillary and eye surgery where the rigidity of airway devices can obscure surgical access. The wire-reinforced tube permits flexion without kinking and can be moved at any time mid-procedure without concern for the cessation of gas flow.



Make an inquiry

Code	Description	Size	Weight	Box Qty.
8005001	Solus Flexible, wire-reinforced laryngeal mask airway	5 Large adult	70+kg (154lbs)	20 
8004001	Solus Flexible, wire-reinforced laryngeal mask airway	4 Medium adult	50–70kg (110–154lbs)	20 
8003001	Solus Flexible, wire-reinforced laryngeal mask airway	3 Small adult	30–50kg (65– 110lbs)	20 
8025001	Solus Flexible, wire-reinforced laryngeal mask airway	2.5 Large pediatric	20–30kg (44–65lbs)	20 
8002001	Solus Flexible, wire-reinforced laryngeal mask airway	2 Small pediatric	10–20kg (22–44lbs)	20 

 Sterile

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 24.7 / IS6.1 US • Issue 15 06.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan  
Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA