


# Solus™ Curve laryngeal mask airway

The laryngeal mask with a preformed airway tube



Airway Management ▪ Airway Devices



## Solus™ Curve laryngeal mask airways

Solus Curve is our laryngeal mask for those who prefer the insertion characteristics of a device with a curved airway tube.

As with all our laryngeal masks, Solus Curve is single use, latex free and supplied sterile.

### The curve

Anatomical curve to aid insertion

### Low friction material

Atraumatic insertion and removal

### Classic cuff shape

Renowned for reliability

### Essential information

Clearly visible for quick reference



### Integral inflation line

Reduces the risk of damage while in use

### High quality valve

To ensure continual cuff integrity

▶ Video available at [us.intersurgical.com](http://us.intersurgical.com)



Code	Description	Size	Weight	Box Qty.
8005050	Solus Curve, large adult, laryngeal mask airway	5	70+kg	20
8004050	Solus Curve, medium adult, laryngeal mask airway	4	50–70kg	20
8003050	Solus Curve, small adult, laryngeal mask airway	3	30–50kg	20

 All supplied sterile

IS6.12 US • Issue 1 02.17