




Solus™ Curve laryngeal mask airway

The laryngeal mask with a preformed airway tube



Airway Management ▪ Airway Devices



Quality, innovation and choice

Interact with us



www.intersurgicalinc.com

Solus™ Curve laryngeal mask airways

Solus Curve is our laryngeal mask for those who prefer the insertion characteristics of a device with a curved airway tube.

As with all our laryngeal masks, Solus Curve is single use, latex free and supplied sterile.

Essential information

Clearly visible for quick reference



The curve

Anatomical curve to aid insertion

Low friction material

Atraumatic insertion and removal

Classic cuff shape


Renowned for reliability

Integral inflation line

Reduces the risk of damage while in use

High quality valve

To ensure continual cuff integrity



Code	Description	Size	Weight	Box Qty.
8005050	Solus Curve, laryngeal mask airway	5 Large adult	70+kg (154lbs)	20 
8004050	Solus Curve, laryngeal mask airway	4 Medium adult	50–70kg (110–154lbs)	20 
8003050	Solus Curve, laryngeal mask airway	3 Small adult	30–50kg (65–110lbs)	20 

[Make an inquiry](#)

[Watch the video](#)

 Sterile

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 181.1 / IS6.12 US • Issue 4 07.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan
Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA