



Intersurgical EcoLite™ adult medium concentration oxygen mask

One of the lightest and most comfortable masks available



Oxygen and Aerosol Therapy • Oxygen Therapy

Intersurgical EcoLite™ adult medium concentration oxygen mask

Welcome to the Comfort Zone

Following the outstanding success of the original Eco adult oxygen masks, with millions sold around the world, we are pleased to announce the arrival of the next generation, Intersurgical EcoLite adult medium concentration oxygen masks.

Incorporating new features such as an extended chin with two integral seals, the new mask ensures a comfortable fit to a wider range of face shapes and sizes.

One of the world's lightest and most comfortable oxygen masks has just become even more comfortable.

Comfortable for the patient

The unique features of the original Eco oxygen mask have been retained and include the following:

Incurved nose seal

Conforms to different nose shapes designed to prevent oxygen entering patient's eyes

No metal nose clip

MRI compatible

Exhalation ports

Designed to minimize the amount of rebreathed carbon dioxide

A choice of "under chin" or "on chin" positions

Provides a better fit on a wider range of patient face shapes

Comfortable for the environment

Various groups have questioned the use of PVC in medical products and it's impact on the environment. Intersurgical's aim is to reduce the environmental impact of its products and processes. The utilization of new materials and the latest manufacturing technology has resulted in the elimination of PVC from the mask shell resulting in a reduced environmental impact.



us.intersurgical.com/info/intersurgicalecolite

Lighter

50% lighter than a PVC mask, more comfortable for the patient

Greener


The elimination of PVC from the masks has dramatically lowered their environmental impact

Soft face seals


Increased patient comfort

Elastic can be positioned under or over the ears

Below ear position eliminates trauma to top of ears



O₂



Make an inquiry

Code	Description	Tube length	Box Qty.
1136015	Intersurgical EcoLite, adult, medium concentration oxygen mask		55
1135015	Intersurgical EcoLite, adult, medium concentration oxygen mask with tube	2.1m	40

Lower environmental impact product

IS10.16 US • Issue 7 02.18