



# Multipac<sup>®</sup> Anesthetic Breathing Systems

A safe and effective one-a-day breathing circuit system that may be used by up to six patients



Anesthesia ▪ Multipac<sup>®</sup> Anesthetic Breathing Systems





## Multipac® Anesthetic Breathing Systems

Cost savings are a major factor in the modern surgical environment. This not only affects products, but also storage requirements and anesthesia set-up time. Intersurgical Incorporated offers a completely integrated anesthesia system to address all of these concerns.

Multipac is a one-a-day breathing circuit that may be used by up to six patients<sup>1</sup>. By selecting one of our eight Multipac configurations, most circuit needs can be met based on your clinical specifications.

By using Multipac systems, cost savings can be realized through the product, the reduction of storage needs and a reduction in waste costs<sup>2</sup>.

Multipac Anesthetic Breathing Circuits include our high-efficiency Filta-Therm™ HMEFs and patient connections. These are packaged within selected circuits and must be exchanged between each patient.

Please see circuit options:



### Multipac anesthetic breathing circuits

Code	Description	Box Qty.
2158009	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite™ size 4 - anesthetic face mask 7294000), one Compact II™ breathing circuit 118 inches long, with parallel Y-piece and end cuffs, 3 liter latex free reservoir bag and 10 foot monitoring line	10
2158010	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 4 - anesthetic face mask 1515000), one Compact II breathing circuit 118 inches long, with parallel Y-piece and end cuffs, 3 liter latex free reservoir bag and 10 foot monitoring line	10
2900077	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 4 - anesthetic face mask 7294000), one UniFlow™ breathing circuit 63 inches long, with end cuffs, 3 liter latex free reservoir bag and 10 foot monitoring line	10
2900078	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 4 - anesthetic face mask 1515000), one UniFlow breathing circuit 63 inches long, with end cuffs, 3 liter latex free reservoir bag and 10 foot monitoring line	10
2902022	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 4 - anesthetic face mask 7294000), one UniFlow breathing circuit 63 inches long with integral monitoring line and end cuffs, 3 liter latex free reservoir bag and 0.9m monitoring line	10
2902023	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 4 - anesthetic face mask 1515000), one UniFlow breathing circuit 63 inches long with integral monitoring line and end cuffs, 3 liter latex free reservoir bag and 0.9m monitoring line	10
2010235	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 4 - anesthetic face mask 7294000), one Flextube™ breathing circuit 63 inches long with parallel Y-piece and end cuffs, 3 liter latex free reservoir bag and 10 foot monitoring line	10
2010261	Multipac Anesthesia Breathing Circuit (2 patient packs include: one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 4 - anesthetic face mask 1515000), one Flextube breathing circuit 63 inches long with parallel Y-piece and end cuffs, 3 liter latex free reservoir bag and 10 foot monitoring line	10

#### References:

- Palmer L., Benbough J., CAMR (Centre for Applied Microbiology & Research) Biosafety Investigation Unit Quality, Safety and Scientific Resources Division Test Report
- Daggan R., Steinberg D., Larjani G., Gratz I., Zafeiridis A., Goldberg M.E.L Intersurgical's Multipac® Breathing System allows reuse of anesthesia breathing circuits (BC) resulting in cost savings and reduced medical waste.

## Multipac® Patient Packs

Additional options are available by selecting one of our eight Patient Packs that will include our high efficiency Filta-Therm™ HMEF, an anesthesia mask elbow and your choice of anesthetic face masks: ClearLite™ or Economy anesthetic face mask. Packs must be exchanged between each patient.



### Multipac Patient Packs

Code	Description	Box Qty.
72936000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 3 - anesthetic face mask 7293000)	20
72946000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 4 - anesthetic face mask 7294000)	20
72956000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 5 - anesthetic face mask 7295000)	20
72966000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one ClearLite size 6 - anesthetic face mask 7296000)	20
15146000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 3 - anesthetic face mask 1514000)	20
15156000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 4 - anesthetic face mask 1515000)	20
15166000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 5 - anesthetic face mask 1516000)	20
15176000	Multipac Patient Pack (Includes one Filta-Therm HMEF 1941000, one anesthesia mask elbow and one Economy size 6 - anesthetic face mask 1517000)	20




## Anesthetic Face Mask Options

### Option 1: ClearLite™ anesthetic face masks

A range of anesthetic face masks which include seven sizes to meet the clinical needs of most patient types and sizes. The range is available in size 0 neonate to size 6 extra large adult. Due to the enhanced anatomical design, the ClearLite mask offers superior seal performance with minimal pressure, reducing the risk of unnecessary discomfort or potential trauma.

The superior seal will also assist the clinician in providing optimum manual ventilation of the patient in both the anesthetic and emergency situation.

ClearLite is both non-PVC and phthalate free which significantly minimizes its impact on the environment. 

#### Lightweight mask

Provides enhanced grip and comfort in use for the clinician

#### Safe and secure connection

This is achieved between the mask and the breathing system in use. The ISO 22mm tapered connection complies with ISO 5356

#### Anatomical shape

Provides a precise fit for different face types and sizes. Minimal pressure required to achieve an effective seal, reducing the risk of inadvertent trauma and nerve damage

#### Flexible non-inflating cuff

Provides a comfortable and effective airtight seal

#### Crystal clear body


Provides excellent visibility for confirmation that the patient is well perfused

#### Color-coding

The comfortable, flexible seal is color coded to allow quick and easy identification of the correct sized mask

#### Reduced mask volume

Reduces potential rebreathing and risks associated with elevated levels of end tidal CO<sub>2</sub>













#### Low environmental impact

Non-PVC and phthalate free

Learn more about ClearLite™ Anesthetic Masks

[www.intersurgical.com/info/ClearLite](http://www.intersurgical.com/info/ClearLite)



Code	Description	Box Qty.
 7296000	ClearLite, anesthetic face mask, size 6, extra large adult, red seal, 22ID	20 
 7295000	ClearLite, anesthetic face mask, size 5, large adult, orange seal, 22ID	35 
 7294000	ClearLite, anesthetic face mask, size 4, adult, green seal, 22ID	35 
 7293000	ClearLite, anesthetic face mask, size 3, small adult, yellow seal, 22ID	35 
 7292000	ClearLite, anesthetic face mask, size 2, pediatric, white seal, 22ID	25 
 7291000	ClearLite, anesthetic face mask, size 1, infant, grey seal, no hook ring, 15OD	30 
 7290000	ClearLite, anesthetic face mask, size 0, neonate, blue seal, no hook ring, 15OD	40 

 Lower environmental impact product

## Anesthetic Face Mask Options

### Option 2: Economy anesthetic face masks

A range of four economy PVC anesthetic face masks with a soft cushion, clear housing and color-coded hook ring. Available in four sizes from small to extra large adult.



#### Clear flexible housing and cushion

Provides great patient visibility

#### Soft cushion

For patient comfort and an air tight fit



#### Color-coded hook ring

For ease of identification, and integration with other Intersurgical airway products

#### Four sizes available

From small to extra large adult

Code	Description	Box Qty.
1517000	Economy, anesthetic face mask, size 6, extra large adult	25
1516000	Economy, anesthetic face mask, size 5, large adult	30
1515000	Economy, anesthetic face mask, size 4, medium adult	40
1514000	Economy, anesthetic face mask, size 3, small adult	50

# Multipac<sup>®</sup>

## anesthetic breathing systems

Multipac<sup>®</sup> is a **safe** and **effective** one-a-day breathing circuit system that may be used by up to six patients



### Keep

Tube, reservoir bag and monitoring line  
for **duration of day**



### Change

Filter, connector and mask  
for **each patient**





**Intersurgical**

6757 Kinne Street  
East Syracuse, NY 13057

T: 800-828-9633  
F: 315-451-3696

[support@intersurgicalinc.com](mailto:support@intersurgicalinc.com)  
[www.intersurgical.com](http://www.intersurgical.com)

