



Intersurgical Guedel Airways

Even Softer by Design

New ISO 5364:2016 Specified Colors

Product Code	Intersurgical Size	ISO Size	Previous Intersurgical Color		NEW ISO Compliant Color	
			Color	Color	Color	Color
1000035	000	3.5	Pink		Light Green	NEW*
1100050	00	5.0	Blue		Blue	
1110055	0	5.5	Grey		Grey	
1111065	1	6.5	White		Brown	NEW**
1111570	1.5	7.0	Yellow		White	NEW***
1112080	2	8.0	Green		Green	
1113090	3	9.0	Orange		Yellow	NEW****
1114100	4	10.0	Red		Red	
1115120	5	12.0	Purple		Purple	

Softer tip

By continuing to prioritize patient safety, Intersurgical has further increased the softness of the tip of all Guedel airways.

Guedel airways that utilize hard and rigid materials at the tip present a significant threat of causing trauma when inserted. By using an even softer material to form the tip of the One-piece Guedel airway, the potential to cause trauma to the patient's airway upon insertion is further reduced.

