

InterGuide Tracheal Tube Introducers (Bougies) and InterForm Intubating Stylets

For use in the Difficult Airway



Introducing InterGuide and InterForm

The InterGuide and InterForm ranges have been designed to complement endotracheal tubes.


Both product ranges are single use, latex free and supplied sterile in easy-to-open individual packs.

Features and Benefits

InterGuide

InterGuide is a flexible tracheal tube introducer commonly known as a “bougie”. It allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations.

Atraumatic coude tip to reduce the potential for patient trauma



Low-friction surface for easy insertion and operation

Available in three sizes: 6FR, 10FR and 15FR

InterForm

InterForm is a malleable endotracheal stylet, allowing the user to form the endotracheal tube into a desired shape to ease insertion and give more control of the airway device.

Aluminium core (MRI safe) allows user to easily form suitable shapes



Low-friction surface for easy insertion and removal

Available in three sizes: 6FR, 10FR and 14FR

Ordering Information

InterGuide



8070010

Code	Size	Diameter (OD)	Length	Box qty.
8070015	15FR	.50cm	70cm	10
8070010	10FR	.33cm	70cm	10
8070006	6FR	.20cm	53cm	10

InterForm



8080010

Code	Size	Diameter (OD)	Length	Box qty.
8080014	14FR	.47cm	34cm	10
8080010	10FR	.33cm	34cm	10
8080006	6FR	.20cm	27.5cm	10

IS6.5 US • Issue 1 12:13



INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

417 Electronics Parkway, Liverpool NY 13088

T: 800-828-9633

support@intersurgicalinc.com

www.intersurgical.com