



Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	QuadraLite masks
Neonatal	Preterm	< 2	0	 Code: 7191000
	Term	2 – 4		
Infant	3 months	6	1	 Code: 7191000
	6 months	8		
	9 months	9		
Pediatric	1 year	10	2	 Code: 7192000
	2 years	12		
	3 years	14		
	4 years	16		
	5 years	18		
Small Adult	6 years	20	3	 Code: 7193000
	7 years	23		
	8 years	26		
	9 years	28		
	10 years	32		
	11 years	35		
	12 years	39		
Medium Adult	13 years	44	4	 Code: 7194000
	14 years	50		
	15 years	54		
	16 years	58		
Large Adult	16+ years	60 – 120	5	 Code: 7194000
	16+ years	> 120	6	 Code: 7194000