




# InterGuide tracheal tube introducers and InterForm intubating stylets

For use in the difficult airway



Airway Management ▪ Airway Accessories



Quality, innovation and choice

Interact with us




[www.intersurgicalinc.com](http://www.intersurgicalinc.com)

## InterGuide tracheal tube introducer

InterGuide is a flexible tracheal tube introducer commonly known as a bougie. It allows positive location of the trachea and subsequent placement of the endotracheal tube in difficult airway situations. InterGuide is single use, latex free and supplied sterile in easy-to-open individual packs.



### Atraumatic coudé tip

To reduce the potential for patient trauma



**Low-friction surface**  
For easy insertion and operation

Make an inquiry


Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8070015	InterGuide tracheal tube introducer bougie	15FR	5.0mm	700mm	10 
8070010	InterGuide tracheal tube introducer bougie	10FR	3.3mm	700mm	10 
8070006	InterGuide tracheal tube introducer bougie	6FR	2.0mm	530mm	10 

## InterForm intubating stylet

InterForm is a malleable stylet, allowing the user to form the endotracheal tube into a suitable shape to ease insertion and give more control of the airway device. InterForm is single use, latex free and supplied sterile in easy-to-open individual packs.




### Aluminium core

Allows user to easily form suitable shapes



**Low-friction surface**  
For easy insertion and removal

Make an inquiry

Code	Description	Size	Outer Diameter (OD)	Length	Box Qty.
8080014	InterForm intubating stylet	14FR	4.7mm	340mm	10 
8080010	InterForm intubating stylet	10FR	3.3mm	340mm	10 
8080006	InterForm intubating stylet	6FR	2.0mm	275mm	10 

 Sterile

**Please think before you print** Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 106.4 / IS6.5 US • Issue 6 06.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK • Australia • Benelux • Brazil • Canada • China • Colombia • Czech Republic • Denmark • France • Germany • Ireland • Italy • Japan  
Lithuania • Malaysia • Mexico • Norway • Philippines • Portugal • Russia • South Africa • Spain • Sweden • Taiwan • Turkey • USA