




## One-piece Guedel airways

Safer by design






## One-piece Guedel airways

### Safer by design

Our Guedel airway combines two materials with very different physical characteristics into one molding. Our one-piece design eliminates the separate bite block found in all conventional Guedel airways. Loose or detached bite blocks can constitute a real threat to patient safety. Our one-piece design overcomes all of the associated risks, as well as improves safety. The one-piece also incorporates a soft tip, reducing the risk of trauma to the patient.

▶ Video available at [www.intersurgical.com](http://www.intersurgical.com)





### ISO 5364:2016 – Color coding guidance

In the most recent version of ISO 5364:2016 Anesthetic and Respiratory Equipment – Oropharyngeal Airways, a new color coding requirement was implemented, aiming to improve patient safety and provide consistency across products from different manufacturers.

In order to comply with the color coding system specified in ISO 5364:2016, it was necessary for Intersurgical to update four sizes of Guedels. The sizes that have changed are as follows:

- ISO size 3.5; from Pink to Light Green
- ISO size 6.5; from White to Brown
- ISO size 7.0; from Yellow to White
- ISO size 9.0; from Orange to Yellow



Each department will need to consider locally how this change will be implemented as it may take time for existing stock to be used.

*The four sizes of Intersurgical Guedel airways that have had their colors changed to conform to the updated standard.*

Product code	ISO size
1000035	3.5
1100050	5.0
1110055	5.5
1111065	6.5
1111570	7.0
1112080	8.0
1113090	9.0
1114100	10.0
1115120	12.0

Previous Intersurgical color	
Pink	
Blue	
Grey	
White	
Yellow	
Green	
Orange	
Red	
Purple	

New ISO compliant color	
Light Green	<b>NEW</b>
Blue	
Grey	
Brown	<b>NEW</b>
White	<b>NEW</b>
Green	
Yellow	<b>NEW</b>
Red	
Purple	












### Softer tip

By continuing to prioritize patient safety, Intersurgical has further increased the softness of the tip of all Guedel airways.

Guedel airways that utilize hard and rigid materials at the tip present a significant threat of causing trauma when inserted. By using an even softer material to form the tip of the One-piece Guedel airway, the potential to cause trauma to the patient's airway upon insertion is further reduced.



Code	Description	Size	ISO	Color	Box Qty.	
1115120	One-piece Guedel airway	5	12.0	Purple	70	
1114100	One-piece Guedel airway	4	10.0	Red	90	
1113090	One-piece Guedel airway	3	9.0	Yellow	90	
1112080	One-piece Guedel airway	2	8.0	Green	100	
1111570	One-piece Guedel airway	1.5	7.0	White	100	
1111065	One-piece Guedel airway	1	6.5	Brown	80*	
1110055	One-piece Guedel airway	0	5.5	Grey	100*	
1100050	One-piece Guedel airway	00	5.0	Blue	100*	
1000035	One-piece Guedel airway	000	3.5	Light Green	100*	

\*Case Qty.

 Lower environmental impact product

US Lit No. 64.3 • Issue date 08.21