




One-piece Guedel airways

Safer by design



Airway Management ▪ Airway Devices



Quality, innovation and choice

Interact with us



www.intersurgicalinc.com



One-piece Guedel airways


Safer by design

By using the latest manufacturing technology Intersurgical has produced the first major breakthrough in Guedel airway design for thirty years.

Combining two materials with very different physical characteristics into one moulding, our One-piece design eliminates the separate bite block found in all conventional Guedel airways.

Loose or detached bite blocks constitute a real threat to patient safety and have been the subject of several hazard notices in the past. Our One-piece design overcomes all of the associated risks.

As well as improved safety, the One-piece Guedel airway also incorporates a soft tip reducing the risk of trauma to the patient. We also offer a sterile option, providing a further choice.





ISO 5364:2016 – Color coding guidance

In the most recent version of ISO 5364:2016 Anesthetic and Respiratory Equipment – Oropharyngeal Airways, a new color coding requirement was implemented, aiming to improve patient safety and provide consistency across products from different manufacturers.

In order to comply with the color coding system specified in ISO 5364:2016, it was necessary for Intersurgical to update four sizes of Guedels. The sizes that have changed are as follows:

- ISO size 9.0; from Orange to Yellow
- ISO size 7.0; from Yellow to White
- ISO size 6.5; from White to Brown
- ISO size 3.5; from Pink to Light Green



ISO Size 9.0 ISO Size 7.0 ISO Size 6.5 ISO Size 3.5

The four sizes of Intersurgical Guedel airways that have had their colors changed to conform to the updated standard.

Each department will need to consider locally how this change will be implemented as it may take time for existing stock to be used.

Product code	ISO size
1115120	12.0
1114100	10.0
1113090	9.0
1112080	8.0
1111570	7.0
1111065	6.5
1110055	5.5
1100050	5.0
1000035	3.5

Previous Intersurgical color	
Purple	
Red	
Orange	
Green	
Yellow	
White	
Grey	
Blue	
Pink	

New ISO compliant color	
Purple	
Red	
Yellow	NEW
Green	
White	NEW
Brown	NEW
Grey	
Blue	
Light Green	NEW

This color change effects all of Intersurgical’s Guedel Airway products, including bulk-packed options.

Softer tip

By continuing to prioritize patient safety, Intersurgical has further increased the softness of the tip of all Guedel airways.

Guedel airways that utilize hard and rigid materials at the tip represent a significant threat of causing trauma when inserted. By using an even softer material to form the tip of the One-piece Guedel airway, the potential to cause trauma to the patient's airway upon insertion is further reduced.



Code	Description	Color	Size	ISO	Box Qty.
1115120	One-piece Guedel airway		5	12.0	70
1114100	One-piece Guedel airway		4	10.0	90
1113090	One-piece Guedel airway		3	9.0	90
1112080	One-piece Guedel airway		2	8.0	100
1111570	One-piece Guedel airway		1.5	7.0	100
1111065	One-piece Guedel airway		1	6.5	20
1110055	One-piece Guedel airway		0	5.5	25
1100050	One-piece Guedel airway		00	5.0	25
1000035	One-piece Guedel airway		000	3.5	25

[Make an inquiry](#)

[Find out more](#)

[Watch the video](#)

Lower environmental impact product

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

U.S. Lit No. 64.5 / IS6.10 US ▪ Issue 7 07.24



Intersurgical Inc, 2055 Executive Dr. Indianapolis, IN 46241, USA

T: 800-828-9633 F: 315-451-3696 support@intersurgicalinc.com www.intersurgicalinc.com

UK ▪ Australia ▪ Benelux ▪ Brazil ▪ Canada ▪ China ▪ Colombia ▪ Czech Republic ▪ Denmark ▪ France ▪ Germany ▪ Ireland ▪ Italy ▪ Japan
Lithuania ▪ Malaysia ▪ Mexico ▪ Norway ▪ Philippines ▪ Portugal ▪ Russia ▪ South Africa ▪ Spain ▪ Sweden ▪ Taiwan ▪ Turkey ▪ USA